

Liquid Acrylamide Solution for Molecular Biology

Poly-acrylamide Gel Electrophoresis solution Protein and Nucleic Acid Separation.

Severn Biotech's high purity liquid Acrylamide Solutions are made from Ultra-Pure Molecular Biology Grade Acrylamide and Bis-acrylamide to exact formulations. They are specific for use in separating peptides, Proteins and Nucleic Acid molecules. Formulations when polymerised give crystal clear electrophoretic gels with reliable reproducable results offering peace of mind separation of biomolecules.

All solutions are manufactured using 18M0hm Ultra-Pure Molecular Biology Grade Water which is LAL tested for endo-toxins and pyrogen free.

The Acryalmide solutions are 0.2micron filtered and DNase RNase and Protease free.

Each batch is MW marker tested to ensure reproducibility.

The solutions are stable and have a long shelf life > 3 years when stored unopened at room tempearture or refrigerated 4°C.

Liquid Acrylamide Solutions available

Product Code	Volume	
30% w/v Acrylamide	Bis-acrylamide Ratio 37.5:1 Protein Gel Formulation	
20-2100-05	500ml	
20-2100-10	1000ml	
30% w/v Acrylamide B	is Acrylamide Ratio 19:1 DNA/ Protein Gel Formulation	
20-2300-05	500ml	
20-2300-10	1000ml	
30% w/v Acrylamide	Bis Acrylamide Ratio 29:1 Protein Gel Formulation	
20-2600-05	500ml	
20-2600-10	1000ml	
40% w/v Acrylamide	Bis Acrylamide Ratio 19:1 DNA Gel Formulation	
20-2400-05	500ml	
20-2400-10	1000ml	
Acrylamide Only Solut	tion - 40% w/v Acrylamide Formulation	
20-2000-05	500ml	
20-2000-10	1000ml	
Bis-acrylamide Only S	Solution - 2% w/v Bis-acrylamide Formulation	
00 0500 05	500ml	
20-2500-05	oodiiit	

Other Acrylamide Stock solutions and formulations are available on request.

Severn Biotech also Manufactures Electrophoresis Buffers used in PAGE.

10X TBE product Code 20-6000-10 1L Pack size 10X TAE Product Code 20-6001-10 1L Pack size 10X Tris Glycine SDS 20-6400-10 1L pack size

Polymerisation Reagents

TEMED Product code 20-3000-01 Pack size 100ml Volume APS Product Code 20-3001-25 Pack size 25g



+44 (0)1562 825286 info@severnbiotech.com severnbiotech.com