

Severn Biotech Ltd.

RNase Removal Surface Cleaning Spray Solution



for molecular biology Laboratories

Pack Size 500ml

This handy 500ml spray solution removes RNase from laboratory surfaces and is effective against enzyme contamination. This will decontaminate bench tops, instruments, pipettes, lab ware, working well on glass or plastic surfaces.

Available as a ready to use solution simply spray on the area contaminated and wipe away from the surface using distilled or RO purified water. Then wipe dry.

Catalogue Numbers 40-1735-05 Volume 500ml

Store & ship at room temperature. Liquid product. For Research Use only

+44 (0)1562 825286 info@severnbiotech.com severnbiotech.com



Unit 2, Park Lane, Kidderminster, Worcestershire DY11 61J Tel: +44 (0)1562 825286 | Fax: +44 (0)1562 825284

Unit 2, Park Lane, Kidderminster, Worcestershire DY11 6TJ

Tel: +44 (0)1562 825286 | Fax: +44 (0)1562 825284 | E-mail: info@severnbiotech.com | www.severnbiotech.com